Contact Us

972-492-0419

0

4100 Fairway Ct. Suite 600 Carrollton, TX 75010



www.realtimelab.com

Organic Acids Profile

CIDS PROFILE

We break down barriers to personal health management and allow people to live better!

realtimelab.com

The Organic Acids Profile is a metabolic snapshot to help you and your patients determine the best course of action in treatment. This urine test looks at metabolic issues, detox impairments, possible dysbiosis, and neurotransmitter metabolites.


Our test measures 36 organic acid markers to create a unique patient biochemistry profile that reflects how well the body is obtaining and utilizing nutrients.

One morning urine sample would translate the body's biochemistry by analyzing the reflection of the absorption of the nutrients gained from your diet and showing any imbalances in vitamins, minerals, and other cofactors that can affect the body's biochemical needs.

Glycolysis


Pyruvate and Lactate are products of the glycolysis cycle (Figure 1) [1] and are key components of many metabolic processes. High values could indicate problems, either inborn or other causes, in multiple different pathways. There are at least five enzymes in this pathway, where a lack of activity could cause elevations of one or more of the markers in this section[1].

Lactate elevation can be caused by exercise [2].However, several toxins can inhibit mitochondrial ATP synthase. This inhibition can lead the body to produce more lactic acid [3]



Fatty Acid Oxidation, Ketone Metabolites, and Protein Metabolism


Succinyl-CoA is the major endpoint of this pathway, which is used in the TCA cycle (see figure 2). Elevations in these markers can be caused by ketogenic diets or exercise. Extreme elevations could be indicative of Maple Syrup urine disease or branched chain keto aciduria [8]. Branched chain amino acids (BCAA) function as metabolic regulators. However, higher levels have been associated with oxidative stress and mitochondrial dysfunction. These issues can lead to insulin resistance in some individuals [9]



Neurotransmitter Metabolism

Neurotransmitters are produced in multiple areas of the body from the central nervous system (CNS), the peripheral nervous system (PNS), and those produced in other parts of the body such as the renal system (see figure 3). The correlation between the CNS neurotransmitter levels and excreted levels are still being researched, although several studies have suggested a link [10, 11]. However, the enzymes that are present in the CNS use the same cofactors. Elevations in precursors could indicate problems. The ratio between vanilmandelate and homovanillate could indicate low expression or activity of dopamine beta hydroxylase or other enzymes in the pathway [12].


figure 3



Detoxification


Glutathione (GSH) is a non-protein tripeptide consisting of three amino acids: glutamate, cysteine, and glycine. (see figure 4) GSH is eventually decomposed into pyroglutamic acid. High amounts of pyroglutamic acid is an indicator of high oxidative stress [13]. Large amounts of pyroglutamic acid could indicate exposure to mycotoxins or other environmental pollutants. Elevated amounts could also indicate a glycine deficiency [14]. Extremely low amounts could indicate a deficiency in the production of glutathione which would cause an individual to become more sensitive to toxins [15].

figure 4



Citric Acid Cycle

Now known as the TCA cycle (see figure 5) but is also called the Krebs' cycle. This cycle is the predominant producer of ATP in most cells[4]. Dietary fruit consumption of orange, grapefruit, and lemon juices can increase urinary citrate levels [5]. Certain levels in this pathway can be increased by either LPS stimulation or toxins [6, 7].



Bacterial Metabolism

Dysbiosis is a common comorbidity for patients suffering from chronic illness. Identifying an overgrowth of pathogenic bacteria can assist in forming a treatment plan. (see figure 6) Markers such as Para-Hydroxyphenlacetate can indicate an overgrowth of Clostridium difficile [16].





figure 5

Sample Report



(uging creatine) Range Perc Ketone Metabolites 20. Alpha-Hydroxybutyrate (H) 4.52 < 1.24 919 21. Beta-Hydroxybutyrate (H) >ULOQ < 8.09 N/A Markers of Neurotransmitter Metabolism 22. Vanilmandelate (H) 4.16 < 3.64 909 23. Homovanillate 3.19 < 6.66 489 24. 5-Hydroxyindoleacetate 2.27 1.17 - 8.06 379 25. Quinolinate 2.74 < 5.37 349 26. Kynurenate 1.25 < 2.49 339 Markers of Detoxification 27.7 1.55 909 28. Orotate 0.42 < 1.04 249 29. Pyroglutamate 34.05 14.58 - 37.47 919 30. Benzoale < LLOQ < 6.87 N/A	3.19
20. Alpha-Hydroxybutyrate (H) 4.52 < 1.24 919 21. Beta-Hydroxybutyrate (H) >ULOQ < 8.09 N/A Markers of Neurotransmitter Metabolism 22. N/A 22. Vanilmandelate (H) 4.16 < 3.64 909 23. Homovanillate 3.19 < 6.66 489 24. 5-Hydroxyindoleacetate 2.27 1.17 - 8.06 379 25. Quinolinate 2.74 < 5.37 349 26. Kynurenate 1.25 < 2.49 339 Markers of Detoxification 27. Para-Hydroxyphenyllactate (H) 1.80 < 1.55 909 28. Orotate 0.42 < 1.04 249 249 399 28. Orotate 0.42 < 1.04 249 30. 919 30. Benzoate < LLOQ < 6.87 N/A	3.19
21. Beta-Hydroxybutyrate (H) >ULOQ < 8.09	3.19
Markers of Neurotransmitter Metabolism 122. Vanilmandelate (H) 4.16 < 3.64	3.19
22. Vanilmandelate (H) 4.16 < 3.64	3.19
23. Homovanillate 3.19 < 6.66	3.19
24. 5-Hydroxyindoleacetate 2.27 1.17 - 8.06 379 25. Quinolinate 2.74 < 5.37	27
25. Quinolinate 2.74 < 5.37	
26. Kynurenate 1.25 < 2.49	274
Markers of Detoxification 27. Para-Hydroxyphenyllactate (H) 1.80 < 1.55	
27. Para-Hydroxyphenyllactate (H) 1.80 < 1.55	125
28. Orotate 0.42 < 1.04 249 29. Pyroglutamate 34.05 14.58 - 37.47 919 30. Benzoate <lloq< td=""> 6.87 N/A</lloq<>	07 -
29. Pyroglutamate 34.05 14.58 - 37.47 919 30. Benzoate <lloq< td=""> < 6.87</lloq<>	
30. Benzoate <lloq 6.87="" <="" a<="" n="" td=""><td>0.42</td></lloq>	0.42
	-
	(NA)-
31. Hippurate 28.66 17.13 - 768.53 109	28.66
Markers of Bacterial Metabolism	
32. Para-Hydroxybenzoate 0.53 < 1.43 609	
33. Para-Hydroxyphenylacetate 15.96 < 26.39 699	-
34. 2-Hydroxyphenylacetate (H) 1.95 < 1.24 93%	-
35. 3-Indoleaœtate 1.71 0.46 - 9.21 45%	\sim
36. Tricarballylate <lloq 1.06="" <="" a<="" n="" td=""><td>171</td></lloq>	171

References

Gray, L.R., S.C. Tompkins, and E.B. Taylor, Regulation of pyruvate metabolism and human disease. Cell Mol Life Sci, 2014. 71(14): p. 2577-604.
Bangsbo, J., F.M. Iaia, and P. Krustrup, Metabolic response and fatigue in soccer. Int J Sports Physiol Perform, 2007. 2(2): p. 111-27.
Ohishi, T., et al., Inhibition of mitochondria ATP synthase suppresses prostate cancer growth through reduced insulin-like growth factor-1 secretion by prostate stromal cells. Int J Cancer, 2020. 146(12): p. 3474-3484.

4. Judge, A. and M.S. Dodd, Metabolism. Essays Biochem, 2020. 64(4): p. 607-647.

5. Barghouthy, Y. and B.K. Somani, Role of Citrus Fruit Juices in Prevention of Kidney Stone Disease (KSD): A Narrative Review. Nutrients, 2021. 13(11).

 Mills, E.L., et al., Succinate Dehydrogenase Supports Metabolic Repurposing of Mitochondria to Drive Inflammatory Macrophages. Cell, 2016. 167(2): p. 457-470 e13.

7. Xie, X., et al., Exposure to polystyrene microplastics causes reproductive toxicity through oxidative stress and activation of the p38 MAPK signaling pathway. Ecotoxicol Environ Saf, 2020. 190: p. 110133.

8.Sgaravatti, A.M., et al., Inhibition of brain energy metabolism by the alpha-keto acids accumulating in maple syrup urine disease. Biochim Biophys Acta, 2003. 1639(3): p. 232-8. 9. Holecek, M., Branched-chain amino acids in health and disease: metabolism, alterations in blood plasma, and as supplements. Nutr Metab (Lond), 2018. 15: p. 33.

10. Lynn-Bullock, C.P., et al., The effect of oral 5-HTP administration on 5-HTP and 5-HT immunoreactivity in monoaminergic brain regions of rats. J Chem Neuroanat, 2004. 27(2): p. 129-38.

11. Chekhonin, V.P., et al., Catecholamines and their metabolites in the brain and urine of rats with experimental Parkinson's disease. Bull Exp Biol Med, 2000. 130(8): p. 805-9.

12. Kundu, P., S. Saha, and G. Gangopadhyay, Kinetics of Allosteric Inhibition of Single Enzyme by Product Molecules. J Phys Chem B, 2020. 124(52); p. 11793-11801.

13. Gamarra, Y., et al., Pyroglutamic acidosis by glutathione regeneration blockage in critical patients with septic shock. Crit Care, 2019. 23(1): p. 162.

14. Jackson, A.A., et al., Urinary excretion of 5-oxoproline (pyroglutamic aciduria) as an index of glycine insufficiency in normal man. Br J Nutr, 1987. 58(2): p. 207-14.

15. Kumar and Bachhawat, Pyroglutamic Acid: throwing light on a lightly studied metabolite. Current science, 2012. 102(2): p. 288-297. 16. Selmer, T., A. Willanzheimer, and M. Hetzel, Propionate CoA-transferase from Clostridium propionicum. Cloning of gene and identification of glutamate 324 at the active site. Eur J Biochem, 2002. 259(1): p. 372-80.